


63-32

CATALOG
AS AD
CTIA
401798



401798



BELL AEROSYSTEMS COMPANY
DIVISION OF BELL AEROSPACE CORPORATION-A  COMPANY

ASTIA
APR 22 1963
RESERVED
ASTIA

BELL AEROSYSTEMS COMPANY
DIVISION OF BELL AEROSPACE CORPORATION

BLR 62-20 (C), Rev. A
April 1963

ANALYTICAL CHEMISTRY

Published and Distributed Under
Contract AF33(657)-8555

BELL AEROSYSTEMS COMPANY
DIVISION OF BELL AEROSPACE CORPORATION

Engineering Laboratories

ANALYTICAL CHEMISTRY

**Bell Laboratory Report
BLR 62-20 (C)
Revision A
April 1963**

Prepared by: Peter Yin
Peter Yin
Laboratory Engineer

Prepared by: William L. Clark
William L. Clark
Unit Leader, Analytical Chemistry

Prepared by: John C. Tynan
John C. Tynan
Group Leader, Analytical Chemistry

Approved by: Henry Ph. Heubusch
Henry Ph. Heubusch
Chief Chemist, Chemistry Section

Approved by: G. F. Kappelt
G. F. Kappelt
Director, Engineering Laboratories

BELL AEROSYSTEMS COMPANY
DIVISION OF BELL AEROSPACE CORPORATION

NOTICE

When US Government drawings, specifications or other data are used for any purpose other than a definitely related government procurement operation, the government thereby incurs no responsibility nor any obligation whatsoever; and the fact that the government may have formulated, furnished, or in any way supplied the said drawings, specifications, or other data is not to be regarded by implication or otherwise, as in any manner licensing the holder or any other person or corporation, or conveying any rights or permission to manufacture, use, or sell any patented invention that may in any way be related thereto. The information contained in this report is for internal use of Bell Aerosystems Company only. It is submitted herewith for information with no guarantee stated or implied as to its accuracy.

Model _____	BELL AEROSYSTEMS COMPANY DIVISION OF BELL AEROSPACE CORPORATION	Page <u>11</u>
Date <u>April 1963</u>		Report <u>BLR 62-20(C)Rev.A</u>

ABSTRACT

Research and development were directed toward the development and improvement of standards or standard procedures for the analyses of rocket propellants. Instrumentation techniques were stressed as a means of obtaining fast, accurate and reliable analyses.

Model _____
Date April 1963

BELL AEROSYSTEMS COMPANY
DIVISION OF BELL AEROSPACE CORPORATION

Page iii
Report BLR 62-20(C), Rev.A

CONTENTS

Section		Page
I	SUMMARY	1
II	INTRODUCTION.....	3
III	RESULTS AND DISCUSSION.....	4
	A. Infrared Methods of Analysis	4
	B. Gas Chromatography.....	6
	C. High Energy Fuel.....	8
	D. Storable Oxidizers	8
	E. Compatibility Studies	10
	F. Compatibility of 100% Dacron.....	10
	G. Analysis of Solid Propellants	10
IV	CONCLUSIONS	15
V	REFERENCES.....	15

Model _____	BELL AEROSYSTEMS COMPANY DIVISION OF BELL AEROSPACE CORPORATION	Page <u>iv</u>
Date <u>April 1963</u>		Report <u>BLR 62-20 (C) Rev. A</u>

ILLUSTRATIONS

Figure		Page
1	Perkin-Elmer Model 421 Grating Infrared Spectrophotometer and its Accessories	5
2	Coleman Carbon - Hydrogen Analyzer	13

TABLES

Number		Page
I	Comparison of Actual Data-H ₂ O Concentrations in 50/50 Fuel Blend as Found by Gas Chromatographic Measurements	7
II	Evaluation of Methods of Analysis for Nitric Oxide in Mixed Oxides of Nitrogen	9
III	Compatibility of Knowlton 071962-A-6 Teflon with Propellants and Flush Fluids	11

Model _____	BELL AEROSYSTEMS COMPANY DIVISION OF BELL AEROSPACE CORPORATION	Page <u>1</u>
Date <u>April 1963</u>		Report <u>BLR 62-20(C), Rev. A</u>

I. SUMMARY

A major highlight of the year has been the acquisition and use of new and advanced instrumentation and accessories. This equipment has greatly increased the capability for the application of modern analytical techniques and has been of value not only in connection with analytical research but also in services rendered to other sections of the Engineering Laboratories Department and to other departments in the Aerospace Rockets and Avionics Division⁸.

Some of the instruments acquired and significant applications are:

- (1) Perkin-Elmer, Model 421, Grating Infrared Spectrophotometer
 - (a) By incorporating a high pressure gas absorption cell in the permeability test apparatus being developed under NASA Contract NAS 7-149 it has been possible to make repeated measurements of the quantity of propellant vapors (MON) permeating through teflon test specimens.
 - (b) Employing specialized sampling accessories, including a microbeam condenser, it was possible to identify the source of contamination which had resulted in the interruption of a gyro test program in the Avionics Division.
 - (c) By making use of an attenuated total reflectance unit it was possible to make a brief preliminary survey of variations in teflon and silicone rubber bladder materials.
 - (d) Using a ten-meter gas absorption cell, a high resolution spectrum of CO₂ gas was recorded for use in checking theoretical calculations in connection with spectroscopic combustion studies.
- (2) Micro Reactor Accessory for the Perkin-Elmer Model 154D Vapor Fractometer.
 - (a) Measurements were made of the desorption of oxygen, nitrogen, carbon dioxide and water from carbon black thermal insulation in an induction furnace. The measurements were made as a function of time and at temperatures up to 380°C.

Model _____	BELL AEROSYSTEMS COMPANY DIVISION OF BELL AEROSPACE CORPORATION	Page <u>2</u>
Date <u>April 1963</u>		Report <u>BLR 62-20(C), Rev.A</u>

(3) Coleman Carbon-Hydrogen Analyzer.

- (a) This instrument has made it possible to analyze synthetic solid propellants within a few hours after preparation. Previously this work was sent to an outside laboratory and at least a week elapsed before the results were available.

Research on the analysis of liquid rocket propellants with particular emphasis on increasing the accuracy and precision of results has continued.

A significant trend during 1962 has been the increasing application of analytical group capability to problems in connection with the purity, contamination, and compatibility of engineering materials. Considerable work has been done with lubricants, hydraulic fluids, solvents, and plastic materials.

Model _____	BELL AEROSYSTEMS COMPANY DIVISION OF BELL AEROSPACE CORPORATION	Page <u>3</u>
Date <u>April 1963</u>		Report <u>BLR 62-20(C), Rev. A</u>

II. INTRODUCTION

Research and development efforts were primarily directed toward the development and/or improvement of standard procedures for the analyses of rocket propellants. Instrumentation techniques were stressed as a means of obtaining fast, accurate and reliable analysis.

Form 0845 Rev. 1248

Model _____	BELL AEROSYSTEMS COMPANY DIVISION OF BELL AEROSPACE CORPORATION	Page <u>4</u>
Date <u>April 1963</u>		Report <u>BLR 62-20(C), Rev. A</u>

III. RESULTS & DISCUSSION

A. INFRARED METHODS OF ANALYSIS

As results of the outstanding performance and versatility of the Perkin-Elmer Model 421 Grating Infrared Spectrophotometer and its accessories (Figure 1), it has contributed strongly to the following achievements of the Analytical Group of the Chemistry Section.

In support of the activities of the Avionic Division, a contamination study on the surface of the slip ring of the gyro assembly was developed and followed up with satisfactory results. Infrared spectroscopic methods have been developed to separate and identify complex organic materials such as epoxy resins, solder flux rosins, and silicone grease. Utilization of the micro sampling accessories available with the instrument provides unlimited possibilities for analytical service to other sections in the future.

Because of the multiplicity of problems involving gaseous and vapor phase samples, no single type of gas cell is adequate for all the measurements. For trace analysis it is necessary to use a long path cell. A 10-meter multiple path cell has been assembled and calibrated.

With the double beam high resolution performance of the Model 421 Infrared Spectrophotometer, a differential technique has been adopted for the identification and semi-quantitative determination of the additive in oil.

A study of the different distillation cuts of the three-component CH_2Cl_2 - H_2O - IPA mixture has been used to evaluate the efficiency of distillation as a method of separation. The composition of each distillation cut was carefully evaluated.

Model _____
Date April 1963

BELL AEROSYSTEMS COMPANY
DIVISION OF BELL AIRSPACE CORPORATION

Page 5
Report BLR 62-20(C), Rev. A




Figure 1. Perkin-Elmer Model 421 Grating Infrared Spectrophotometer and its Accessories

Model _____	BELL AEROSYSTEMS COMPANY DIVISION OF BELL AEROSPACE CORPORATION	Page <u>6</u>
Date <u>April 1963</u>		Report <u>BLR 62-20(C), Rev. A</u>

A new method for the production of intense and sensitive reflection spectra of strongly absorbing substances has been employed. The spectrum thus obtained shows a high degree of similarity to ordinary absorption spectra. This Attenuated Total Reflection (ATR) technique will be useful in many cases in which normal infrared techniques fail due to the difficulties in sample preparation. The special features of this ATR technique find important applications in the solid propellant project where the samples are in a spongy solid state and no solvents are available.

B. GAS CHROMATOGRAPHY

The micro reactor accessory of the Model 154D Vapor Fractometer was installed and operated satisfactorily. The techniques and applications of the micro reactor were critically evaluated. The study of the desorption of H₂O from lamp black as a function of temperature and time has proved the micro reactor to be highly useful for such types of analysis.

An initial attempt to trap the gas chromatographic effluents of the lamp black at various temperatures for further identification by infrared spectroscopy was satisfactory. More study has to be done along this line in order to increase the usage of the micro reactor in the future.

The printing integrator for the Model 154D Vapor Fractometer was installed and checked satisfactorily. A statistical evaluation program was initiated to compare the accuracy of methods of calculation by peak height, printing integrator, and planimeter. The resultant data was programmed in an IBM 704 computer and it was found that the corresponding deviation of each method was within the experimental error (see Table I).

Theoretically, the printing integrator should give the highest accuracy. After an adjustment had been made by the vendor another check was made. H₂O in the range of 0.01 to 0.04 was measured with an average deviation of ± 0.0001 cc.

Model _____
Date April 1963

BELL AEROSYSTEMS COMPANY
DIVISION OF BELL AEROSPACE CORPORATION

Page 7
Report BLR 62-20(C), Rev. A

TABLE I
COMPARISON OF ACTUAL DATA-H₂O CONCENTRATIONS
IN 50/50 FUEL BLEND AS FOUND BY GAS
CHROMATOGRAPHIC MEASUREMENTS

<u>Actual H₂O Present %W</u>	<u>Integrator</u>	<u>Planimeter</u>	<u>Peak Height</u>
0.48	0.52	0.51	0.49
	0.58	0.51	0.49
		0.52	0.50
1.01	0.93	0.99	1.02
	0.92	0.92	0.99
		0.93	0.99
1.61	1.62	1.70	1.64
	1.69	1.58	1.60
2.17	2.13	2.20	2.19
	2.11	2.14	2.16
		2.18	2.15
2.87	2.93	2.87	2.85
	2.91	2.89	2.89
		2.89	2.89
Avg. Deviation	±0.06	±0.04	±0.02

Model _____	BELL AEROSYSTEMS COMPANY DIVISION OF BELL AEROSPACE CORPORATION	Page <u>8</u>
Date <u>April 1963</u>		Report <u>BLR 62-20(C), Rev.A</u>

The molecular sieve V-A column for gas chromatography was found very satisfactory in separating the CO_2 - CO , O_2 and N_2 quantitatively.

C. HIGH ENERGY FUEL

50/50 UDMH and N_2H_4 Fuel Blend - A gas chromatographic field test method of water content in the 50/50 fuel blend was developed as a laboratory procedure for routine use. The direct determination of water content is perhaps the most difficult problem presented, however, the precision and accuracy of this method was found to be ± 0.02 percent.

A near infrared spectrophotometric method has been adopted for the routine determination of the 50/50 UDMH and N_2H_4 content of the fuel blend. The Beckman DK-1 ratio recording spectrophotometer was employed for this work.

D. STORABLE OXIDIZERS

Allied Chemical Company's procedure for the determination of nitric oxide in MON was evaluated (Reference 1). The principle of this method was a direct oxidation of the NO in a 2 inch diameter stainless steel bomb to form N_2O_4 . By the difference of the weights of the stainless steel bomb before and after the oxidation, the percent weight of NO can be calculated.

The precision of the existing melting point and titration methods and the bomb method were evaluated statistically (see Table II).

As indicated in the table, the bomb method has better precision and gives slightly higher results than the other two at the 95 percent confidence level. A routine procedure has been established in the Propellants Laboratory.

To obtain reliable analytical data for the evaluation of a permeability test apparatus, a pair of high pressure infrared gas analysis cells were procured. These cells were modified in order to improve their compatibility with rocket propellant vapors. Calibration runs were made with a series of nitrogen tetroxide/nitrogen gas mixtures up to 200 psig total pressure.

Model

Date

April 1963

BELL AEROSYSTEMS COMPANY

DIVISION OF BELL AEROSPACE CORPORATION

Page

9

Report

BLR 62-20(C), Rev. A

TABLE II

EVALUATION OF NITRIC OXIDE DETERMINATION METHOD

No. of Determinations	Melting Point	Neutralization	Bomb	Bomb
n = 1	9.5	9.72	10.37	-
2	9.8	10.10	10.37	-
3	9.5	9.91	10.59	10.59
4	9.5	10.00	10.51	10.51
5	9.8	10.07	10.54	10.54
6	9.8	10.01	10.49	10.49
7	10.0	10.18	10.52	10.52
8	9.4	10.10	10.51	10.51
9	9.8	10.07	10.51	10.51
10	9.8	10.15	10.47	10.47
11	-	-	10.53	
12	-	-	10.59	
13	-	-	10.56	
\bar{V}	0.039	0.018	0.005	0.0011
S	0.196	0.134	0.073	0.033
95% CL $P = \frac{TS}{\sqrt{n}}$	±.14	±.09	±0.043	±0.042
X	9.69	10.03	10.505	10.517

Model _____	BELL AEROSYSTEMS COMPANY DIVISION OF BELL AEROSPACE CORPORATION	Page <u>10</u>
Date <u>April 1963</u>		Report <u>BLR 62-20(C), Rev. A</u>

Data was obtained which indicated that inert pressurizing gases up to 200 psig do not significantly effect the intensity of the N_2O_4 absorption bands used for quantitative analyses.

E. COMPATIBILITY STUDIES

A study of the effect of curing on the compatibility of Knowlton 071962-A-6 teflon with various oxidizers, fuels, and solvents was undertaken. The specified curing cycle was 5 minutes at 500°F followed by 15 minutes at 600°F.

A preliminary test showed that while the teflon withstood the 500°F temperature, it was badly charred at 600°F. Therefore, the compatibility tests were run after the first curing stage. The results are given in Table III.

F. COMPATIBILITY OF 100% DACRON

The results of a seven-day compatibility test of 100% Dacron in different fluids can be summarized as follows.

The 100% Dacron disintegrated immediately on contact with oxidizers such as IRFNA, MON, N_2O_4 . In the fuels, N_2H_4 and the blend, it broke up within 24 hours. No attack was visibly noticed when the 100% Dacron was contacted with the flush fluids such as methanol, isopropanol, methylene chloride, and inhibited water.

G. ANALYSIS OF SOLID PROPELLANTS

An automatic carbon-hydrogen analyzer (Figure 2) was purchased and checked out for performance, and is being used for the micro determination of carbon and hydrogen contents of solid propellants.

A remote control Schöniger Model 141R apparatus for the determination of boron was also installed and calibrated for laboratory use.

Form 62-20 Rev. 12-62

Model _____		BELL AEROSYSTEMS COMPANY		Page <u>11</u>
Date <u>April 1963</u>		DIVISION OF BELL AEROSPACE CORPORATION		Report <u>BLR 62-20(C), Rev. A</u>
TABLE III				
COMPATIBILITY TESTING				
Fluids	Physical Appearance of Teflon	Physical Appearance of Fluid	Wt. Loss .gm	
1. N ₂ H ₄	Noncured	Light brown-purple (0.4728 gm Sample wt.)	N ₂ H ₄ slightly tinted- yellowish - no residue	-0.0079
	Cured	No discoloration (0.5832 gm)	N ₂ H ₄ discolored yellowish brown - no residue	-0.0133
2. UDMH	Noncured	Light brownish - purple (0.4132 gm)	UDMH clear	-0.0025
	Cured	No discoloration (.2273 gm)	UDMH clear	-0.0033
3. 50/50 Blend	Noncured	Light brown-purple (0.5171 gm)	Clear liquid	+0.0048
	Cured	(0.2756 gm)	Blend discolored only slightly	-0.0040
4. Methanol	Noncured	Light brown - purple (0.4744 gm)	Clear liquid	-0.0100
	Cured	No discoloration (0.2613 gm)	Clear liquid	-0.0043

Model _____		BELL AEROSYSTEMS COMPANY DIVISION OF BELL AEROSPACE CORPORATION		Page <u>12</u>
Date <u>April 1963</u>				Report <u>BLR 62-20(C), Rev. A</u>
TABLE III (CONT)				
Fluids	Physical Appearance of Teflon		Physical Appearance of Fluid	Wt. Loss gm
5. Isopropanol	Noncured	Light brown - purple (0.4398 gm)	Clear liquid	-0.0086
	Cured	No discoloration (0.5113 gm)	Clear liquid	-0.0043
6. Methylene Chloride	Noncured	Light brown - purple (0.4063)	Clear liquid	-0.0107
	Cured	Clear (0.6082)	Clear liquid	-0.0017
7. Chromic Acid	Noncured	Light brown - purple (0.4744 gm)	Clear liquid	-0.0034
	Cured	No discoloration (0.2706)	Clear liquid	-0.0054
8. N ₂ O ₄ (Cured Only)		Bleach coupons white (0.6115 gm)	Liquid evaporated within 24 hours leaving a thin layer of residue at bottom of beaker	0.07
		Bleached coupons white (0.5509 gm)	Liquid evaporated within 24 hours heaving a thin layer of residue at the bottom of beaker	0.0583
9. MON (Cured Only)				
10. IRFNA (Cured Only)		Teflon appeared to shred and loosened binder immediately	No visible change	-

Model _____
Date April 1963

BELL AEROSYSTEMS COMPANY
DIVISION OF BELL AEROSPACE CORPORATION

Page 13
Report BLR 62-20(C), Rev. A




Figure 2. Coleman Carbon-Hydrogen Analyzer

Model _____	BELL AEROSYSTEMS COMPANY DIVISION OF BELL AIRSPACE CORPORATION	Page <u>14</u>
Date <u>April 1963</u>		Report <u>BLR 62-20(C), Rev. A</u>

By taking advantage of the variety of operating modes of the 421 Spectrophotometer it has been possible to obtain high resolution spectra of PEH. This has not previously been possible. It is believed that these spectra will make it possible to make a more complete correlation between the infrared absorption and the polymer structure.

Model _____	BELL AUROSYSTEMS COMPANY DIVISION OF BELL AIRSPACE CORPORATION	Page <u>15</u>
Date <u>April 1963</u>		Report <u>BLR 62-20(C), Rev. A</u>

IV. CONCLUSIONS

The availability of new, modern instrumentation and laboratory equipment has made it possible for the Analytical Group to accomplish its goal of increased accuracy, precision, and rapidity of analytical procedures. In addition, the scope of applications has been broadened to the point that the analysis of materials other than propellants has become an important phase of the operations with no significant decrease in propellant chemistry.

V. REFERENCES

1. Analytical Procedures for Mixed Oxides of Nitrogen, Allied Chemical Corporation, Hopewell, Virginia.